

Supris[®]

Retropubic Sling System

Purposefully designed

Female Stress Urinary Incontinence





Retropubic Sling System

Coloplast partnered with physicians to design the patented Supris® Retropubic Sling System which consists of the Supris implantable midurethral sling and disposable introducers that allow for either **a retropubic top-down or bottom-up approach**. It is indicated for the surgical treatment of female stress urinary incontinence (SUI) resulting from urethral hypermobility and/or intrinsic sphincter deficiency.¹

Dual-use
introducers

Distinct curve for use in both
a retropubic top-down or
bottom-up approach

Predictable
tension

Mesh facilitates
positioning
during surgery¹

Macroporous
design

allows for optimal
tissue integration¹

Coloplast SUI Mesh

Predictable tension mesh that facilitates positioning during surgery. Macroporous design allows for optimal tissue integration and provides a surface area as a backboard to support the urethra.¹ It may look a little different but the reasons are purposefully clear:

Pore size
About 5x the minimum

Thickness
About half as thin as others on
the market: ≤ .30 mm

Density
Lower than others
on the market: 80 g/m²

Elasticity
Low: 7.5%

Fiber size
Lower than others
on the market: .08 mm

Total material/mass
Low amount of total
material implanted

To Order Call Toll Free 800.258.3476

This product may be ordered directly from Coloplast.

Supris® Retropubic Sling System	Purchase UOM	Qty per UOM	Item
60 cm Sling	EA	1	93-4450
Retropubic Introducers (top-down or bottom-up)	EA	2	

Coloplast Delivers Confidence



Scan to learn more
about **Supris®**

It all started in 1954 with Elise Sorensen, a nurse who wanted to help her sister regain confidence and control. After having an ostomy operation, her sister was afraid to go out in public for fear of leakage. Elise worked with engineer Aage Louis-Hansen to develop the world's first adhesive ostomy bag for greater protection and control. Today, Coloplast is an environmentally conscious global leader in ostomy care, wound and skin care, and urology care products that improve lives for millions of people around the world.

SUPRIS® RETROPUBIC SLING SYSTEM BRIEF STATEMENT

Indications:

The Supris Retropubic Kit consists of the Supris implantable midurethral support sling and disposable introducers for placement using a "top-down" or "bottom-up" retropubic surgical approach. The Supris sling and introducers are indicated for the surgical treatment of female stress urinary incontinence (SUI), resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

Contraindications:

It is the responsibility of the physician to advise the prospective patients or their representatives, prior to surgery, of the contraindications associated with the use of this product. The Supris Retropubic Kit is contraindicated for use in patients with the following conditions:

- Pregnancy or desire for future pregnancy
- Potential for further growth (e.g., adolescents)
- Known active urinary tract infection and/or infection in operative field
- Taking anti-coagulant therapy
- Abnormal urethra (e.g., fistula, diverticulum)
- Intraoperative urethral injury
- Any condition, including known or suspected pelvic pathology, which could compromise implant or implant placement
- Sensitivity/allergy to polypropylene

Warnings and Precautions:

It is the responsibility of the physician to advise the prospective patients or their representatives, prior to surgery, of the warnings and precautions associated with the use of this product and the associated surgical risks.

Warnings:

The Supris Retropubic Kit should only be used by physicians familiar with the surgical procedures and techniques involving transvaginal placement of non-absorbable, synthetic mesh slings and who have adequate education and experience in the treatment of female SUI.

A thorough assessment of each patient should be made to determine the suitability of a synthetic mesh sling procedure.

The patient should be counseled that alternative incontinence treatments may be appropriate, and the reason for choosing a mesh sling procedure should be explained.

Obtain patient consent prior to surgery and ensure that the patient has an understanding of the postoperative risks and potential complications of transvaginal mesh sling surgery.

Patient counseling should include a discussion that the sling to be implanted is a permanent implant and that some complications associated with the implanted mesh sling may require additional surgery; repeat surgery may not resolve these complications. Serious adverse tissue responses or infection may require removal of mesh, and complete removal of the sling may not always be possible. Individuals may have varying degrees of collagen laydown that may result in scarring.

As with all surgical procedures, patients with certain underlying conditions may be more susceptible to postoperative bleeding, impaired blood supply, compromised/delayed healing, or other complications and adverse events.

The risks and benefits of using Supris should be considered in patients.

Any future pregnancy could negate the benefits of this surgical procedure. Patients should report any bleeding, pain, abnormal vaginal discharge or sign of infection that occur at any time.

Do not use product that has damaged or opened packaging, or has expired, as sterility may be compromised.

The procedure to insert the Supris sling requires good knowledge of pelvic anatomy and the correct use of the introducer needles in order to avoid damage to adjacent anatomical structures.

Potential Complications

Adverse events are known to occur with transvaginal synthetic sling procedures and implants. Adverse events following mesh implantation may be de novo, persistent, worsening, transient, or permanent.

Adverse events may include but are not limited to: abscess (acute or delayed), adhesion/scar formation, allergy, hypersensitivity or other immune reaction, bleeding, hemorrhage or hematoma, bowel obstruction, dehiscence,

delayed wound healing, extrusion, erosion or exposure of mesh sling into the vagina or other structures or organs, fistula formation, infection, inflammation (acute or chronic), local irritation, necrosis, de novo and/or worsening dyspareunia, neuromuscular symptoms (acute or chronic), pain (acute or chronic), partner pain and/or discomfort during intercourse, perforation or injury of soft tissue (e.g., muscles, nerves, vessels), structures, or organs (e.g., bone, bladder, urethra, ureters, vagina), seroma, sling migration, bladder storage dysfunction (e.g., increased daytime frequency, urgency, nocturia, overactive bladder, urinary incontinence), ureteral obstruction, urinary tract infection, voiding symptoms (e.g., dysuria, urinary retention, incomplete emptying, straining, positional voiding, weak stream), granulation tissue formation, palpable mesh (patient and/or partner), sexual dysfunction, vaginal discharge (abnormal) and vaginal scarring or tightening.

The occurrence of these events may require one or more revision surgeries, including removal of the sling.

Complete removal of the sling may not always be possible, and additional surgeries may not always fully correct the complications.

There may be unresolved pain with or without mesh sling explantation.

The information provided is not comprehensive with regard to product risks. For a comprehensive listing of indications, contraindications, warnings, precautions, and adverse events refer to the product's Instructions for Use. Alternatively, you may contact a Coloplast representative at 1-800-258-3476 and/or visit the company website at www.coloplast.com.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

04.2021
PM-04431

References

1. Data on file at Coloplast.