

Aris[®]

Transobturator Sling System

Purposefully designed

Female Stress Urinary Incontinence




Aris[®]

Transobturator Sling System

Coloplast partnered with physicians to design the patented Aris[®] Transobturator Sling. Each system includes one flat curve and one set of helical introducers, providing physicians the option to choose their preferred transobturator (TO) approach for each patient.





Low complications

4% Chronic groin pain, with 1.4% requiring medical or surgical treatment; 8% continence reoperation rate at 9-year follow-up¹

Coloplast SUI Mesh

Predictable tension mesh that facilitates positioning during surgery. Macroporous design allows for optimal tissue integration and provides a surface area as a backboard to support the urethra.⁵ It may look a little different but the reasons are purposefully clear:

Pore size

About 5x the minimum

Thinness

About half as thin as others on the market: $\leq .30$ mm

Density

Lower than others on the market: 80 g/m^2

Elasticity

Low: 7.5%

Fiber size

Lower than others on the market: .08 mm

Total material/mass

Low amount of total material implanted

To Order Call Toll Free 800.258.3476

This product may be ordered directly from Coloplast.

Aris® Transobturator Sling System	Purchase UOM	Qty per UOM	Item
60 cm Sling	EA	1	93-4400
Flat Inducer	EA	1	
Helical Inducers	EA	2	



Coloplast Delivers Confidence

It all started in 1954 with Elise Sorensen, a nurse who wanted to help her sister regain confidence and control. After having an ostomy operation, her sister was afraid to go out in public for fear of leakage. Elise worked with engineer Aage Louis-Hansen to develop the world's first adhesive ostomy bag for greater protection and control. Today, Coloplast is an environmentally conscious global leader in ostomy care, wound and skin care, and urology care products that improve lives for millions of people around the world.

ARIS® TRANSOBTURATOR KIT BRIEF STATEMENT

Indications:

The Aris Transobturator Kit consists of the Aris implantable midurethral support sling and disposable introducers. The Aris sling and introducers are indicated for the surgical treatment of all types of stress urinary incontinence (SUI) and for female urinary incontinence resulting from urethral hypermobility and/or intrinsic sphincter deficiency. The introducer is a surgical instrument designed to assist in correct placement of a sling.

Contraindications:

It is the responsibility of the physician to advise the prospective patients or their representatives, prior to surgery, of the contraindications associated with the use of this product. The Aris Transobturator Kit is contraindicated for use in patients with the following conditions:

- Pregnancy or desire for future pregnancy
- Potential for further growth (e.g., adolescents)
- Known active urinary tract infection and/or infection in operative field
- Taking anti-coagulant therapy
- Abnormal urethra (e.g., fistula, diverticulum)
- Intraoperative urethral injury
- Any condition, including known or suspected pelvic pathology, which could compromise implant or implant placement
- Sensitivity/allergy to polypropylene

Warnings and Precautions:

It is the responsibility of the physician to advise the prospective patients or their representatives, prior to surgery, of the warnings and precautions associated with the use of this product and the associated surgical risks.

Warnings:

The Aris Transobturator Kit should only be used by physicians familiar with the surgical procedures and techniques involving transvaginal placement of non-absorbable, synthetic mesh slings and who have adequate education and experience in the treatment of female SUI.

A thorough assessment of each patient should be made to determine the suitability of a synthetic mesh sling procedure.

The patient should be counseled that alternative incontinence treatments may be appropriate, and the reason for choosing a mesh sling procedure should be explained.

Obtain patient consent prior to surgery and ensure that the patient has an understanding of the postoperative risks and potential complications of transvaginal mesh sling surgery.

Patient counseling should include a discussion that the sling to be implanted is a permanent implant and that some complications associated with the implanted mesh sling may require additional surgery; repeat surgery may not resolve these complications. Serious adverse tissue responses or infection may require removal of mesh, and complete removal of the sling may not always be possible. Individuals may have varying degrees of collagen laydown that may result in scarring.

As with all surgical procedures, patients with certain underlying conditions may be more susceptible to postoperative bleeding, impaired blood supply, compromised/delayed healing, or other complications and adverse events.

The risks and benefits of using Aris should be considered in patients.

Any future pregnancy could negate the benefits of this surgical procedure. Patients should report any bleeding, pain, abnormal vaginal discharge or sign of infection that occur at any time.

Do not use product that has damaged or opened packaging, or has expired, as sterility may be compromised.

The procedure to insert the Aris sling requires good knowledge of pelvic anatomy and the correct use of the introducer needles in order to avoid damage to adjacent anatomical structures.

Potential Complications

Adverse events are known to occur with transvaginal synthetic sling procedures and implants. Adverse events following mesh implantation may be de novo, persistent, worsening, transient, or permanent.

Adverse events may include but are not limited to: abscess (acute or delayed), adhesion/scar formation, allergy, hypersensitivity or other immune reaction, bleeding, hemorrhage or hematoma, dehiscence, delayed wound healing, extrusion, erosion or exposure of mesh sling into the vagina or other structures or organs, fistula formation, infection, inflammation (acute or chronic), local irritation, necrosis, de novo and/or worsening dyspareunia, neuromuscular symptoms (acute or chronic), pain (acute or chronic), partner pain and/or discomfort during intercourse, perforation or injury of soft tissue (e.g., muscles, nerves, vessels), structures, or organs (e.g., bone, bladder, urethra, ureters, vagina), seroma, sling migration, bladder storage dysfunction (e.g., increased daytime frequency, urgency, nocturia, overactive bladder, urinary incontinence), ureteral obstruction, urinary tract infection, voiding symptoms (e.g., dysuria, urinary retention, incomplete emptying, straining, positional voiding, weak stream), granulation tissue formation, palpable mesh (patient and/or partner), sexual dysfunction, vaginal discharge (abnormal) and vaginal scarring or tightening.

The occurrence of these events may require one or more revision surgeries, including removal of the sling.

Complete removal of the sling may not always be possible, and additional surgeries may not always fully correct the complications.

There may be unresolved pain with or without mesh sling explantation.

The information provided is not comprehensive with regard to product risks. For a comprehensive listing of indications, contraindications, warnings, precautions, and adverse events refer to the product's Instructions for Use. Alternatively, you may contact a Coloplast representative at 1-800-258-3476 and/or visit the company website at www.coloplast.com.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

PM-04430 09/2021

References

1. Karmakar D, Mostafa A, Abdel-Fattah M. Long-term outcomes of transobturator tapes in women with stress urinary incontinence: E-TOT randomised controlled trial. *BJOG* 2017;124:973-981.
2. Nilsson CG, Palva K, Aarnio R, Morcos E, Falconer C. Seventeen years' follow-up of the tension-free vaginal tape procedure for female stress urinary incontinence. *Int Urogynecol J* 2013;24: 1265-9.
3. Ulrich D, Tammaa A, Holbfer S, Trutnovsky G, Bjelic-Radiscic V, Tamussino K, et al. Ten-year followup after tension-free vaginal tape-obturator procedure for stress urinary incontinence. *J Urol* 2016;196:1201-6.
4. Serati M, Braga A, Athanasiou S, Tommaselli GA, Caccia G, Torella M, et al. Tension-free vaginal tape-obturator for treatment of pure urodynamic stress urinary incontinence: efficacy and adverse effects at 10-year follow up. *Eur Urol* 2016; DOI: 10.1016/j.eururo.2016.
5. Data on file.