

Ureteral Dilators

Coloplast Ureteral dilators were designed to facilitate ureteral dilatation: they are radiopaque for monitoring the advancement under fluoroscopy, and their diameter is indicated at each tip for convenience.

- Length: 48 cm
- Compatible with guidewire up to 0.038"
- Diameter indicated at each tip
- Packaged in box of 10

Diameter (Ch)	Diameter (Fr)	Length (cm)	Item
8	10	48	RBD010
12	14	48	RBD014



Ureteral Dilator Brief Statement

Indications

Dilatation of the ureter during ureteroscopy procedures.

Contraindications

Same as for ureteral dilatation.

Any contraindication to endoscopy.

Any known allergies to the medical device materials.

Warnings and Precautions

This type of device must only be used by trained and experienced professionals.

Reuse of this single use product may create a potential risk to the user.

Potential Complications

Hematuria due to local ureteral bleeding may occur.

As with any dilation procedure, the use of a ureteral dilator may cause traumatic tissue lesions.

There is a risk of injury to the urinary tract, particularly if the instructions for use have not been complied with, and especially if the dilator has been positioned without radioscopic control.

The information provided is not comprehensive with regard to product risks. For a comprehensive listing of indications, contraindications, warnings, precautions, and adverse events refer to the product's Instructions for Use. Alternatively, you may contact a Coloplast representative at 1-800-258-3476 and/or visit the company website at www.coloplast.com.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Minneapolis, MN

12/7/21

PM-19627