

Torosa[®]

Saline-Filled Testicular Prosthesis

The only **FDA-approved** testicular implant

Torosa[®] is the only FDA-Approved testicular implant available to aid in creating and restoring a more natural look and feel of the testicles in the scrotum.


Torosa is available in multiple sizes to meet the needs of adult, adolescent and pediatric patients.

Our single-use sizing guide will assist you in determining the best size for your patient.




Torosa[®]


Saline-Filled Testicular Prosthesis




Extra Small: 7 cc



Small: 11 cc



Medium: 16 cc



Large: 20 cc

Ordering Information

Coloplast Interventional Urology Surgical Support **800-258-3476**

Order Number	Size	Width	Height	Volume
450-1323	Extra Small	2.2 cm	3.0 cm	7 cc
450-1325	Small	2.5 cm	3.5 cm	11 cc
450-1327	Medium	2.7 cm	4.0 cm	16 cc
450-1329	Large	2.9 cm	4.5 cm	20 cc

For additional Torosa resources including educational videos, filling instructions, coding guide and additional ordering information, visit ColoplastMD.com or [scan the QR Code](#).



TOROSA® SALINE-FILLED TESTICULAR PROSTHESIS BRIEF STATEMENT

Indications: The Coloplast TOROSA Saline-Filled Testicular Prosthesis is intended for use when cosmetic testicular replacement is indicated i.e., in the case of agenesis or following the surgical removal of a testicle.

Contraindications: The implantation of testicular prostheses is contraindicated in the presence of infection or untreated neoplasm.

Warnings: This device contains solid silicone elastomer. The risks and benefits of implanting this device in patients with lupus (e.g., SLE or DLE), scleroderma (e.g., progressive systemic sclerosis), myasthenia gravis, or documented sensitivity to silicone should be carefully considered. The issue of the possible relationship between silicone and various diseases has been and continues to be the subject of scientific and medical debate. Sepsis or hemorrhage may result from the placement of any foreign object in the body.

Excessive fibrous capsular formation or contracture may occur around any implant placed in contact with soft tissues.

Precautions: Each prosthesis should be checked for patency prior to surgery and continuously monitored throughout the surgical procedure to ensure that the structural integrity of the implant is not compromised in any way. The action of drugs (such as antimicrobials, chemotherapy agents or steroids) in contact with the prosthesis has not been tested by the manufacturer, and their use cannot be recommended.

Each physician who chooses to use drugs in combination with this prosthesis must assure compatibility of the drug with silicone elastomer.

A thorough preoperative consultation should include a discussion between the patient and physician of all available treatment options and their risks and benefits.

Implantation of the TOROSA Saline-Filled Testicular Prosthesis may be difficult or impossible in patients with inadequate scrotal tissue to cover the prosthesis, patients who have undergone prior pelvic radiation therapy, or patients whose wound healing abilities are compromised (e.g., uncontrolled diabetes, poor circulation).

Potential Complications: Potential complications include pain and discomfort.

The information provided is not comprehensive with regard to product risks. For a comprehensive listing of indications, contraindications, warnings, precautions, and adverse events refer to this product's Instructions for Use.

Alternatively, you may contact a Coloplast Representative at 1-800-258-3476 and/or visit the company website at www.coloplast.us.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

[Ostomy Care / Continance Care / Wound & Skin Care / Interventional Urology](#)


Coloplast Corp. Minneapolis, MN 55411 / Interventional Urology Surgical Support 1-800-258-3476

www.coloplast.us The Coloplast logo is a registered trademark of Coloplast A/S.


© 2021 Coloplast Corp. All rights reserved.



PM-00124 07.21



Large: 20 cc



Medium: 16 cc



Small: 11 cc



Extra Small: 7 cc

Torosa®
Saline-Filled Testicular Prosthesis