

Vortek[®] Hydro

Hydrophilic-Coated PEBA Ureteral Stent

*Your new go-to for
performance and value*

Dual Layer Design

- Stiff two-layer structure aids in placement and resists compression
- Thermosensitive outer layer softens at body temperature
- Rigid internal polymer allows for graceful glide over guidewire

Slick Hydrophilic-Coating

- Facilitates insertion and advancement
- Improves glide between the stent and epithelial layer of ureter
- Double layer structure allows for easy insertion over guidewire

Efficient Urinary Drainage

- Drains urinary tract over ureteral obstacles
- Drainage capacity of 4.8 Fr is similar to 6 Fr
- Radiopaque & indicated for 6 month indwell

For ordering and additional information please contact your sales representative or Coloplast Interventional Urology Surgical Support at [800-258-3476](tel:800-258-3476)



Ordering Information

Coloplast Interventional Urology Surgical Support 800-258-3476

French	Length	Order Number Without Guidewire	Order Number Orchestra™ Guidewire	Order Number Seldinger Guidewire
4.8	22 cm	BCFD47	BCFS47	BCFA47
	24 cm	BCFD43	BCFS43	BCFA43
	26 cm	BCFD44	BCFS44	BCFA44
	28 cm	BCFD45	BCFS45	BCFA45
6	22 cm	BCFD67	BCFS67	BCFA67
	24 cm	BCFD63	BCFS63	BCFA63
	26 cm	BCFD64	BCFS64	BCFA64
	28 cm	BCFD65	BCFS65	BCFA65
	30 cm	BCFD66	BCFS66	BCFA66
7	22 cm	BCFD77	BCFS77	BCFA77
	24 cm	BCFD73	BCFS73	BCFA73
	26 cm	BCFD74	BCFS74	BCFA74
	28 cm	BCFD75	BCFS75	BCFA75
	30 cm	BCFD76	BCFS76	BCFA76

Indications for Use: Drainage of the upper urinary tract over fistulas or ureteral obstacles. Healing of the ureter.

Warnings: Reuse of this single use product may create a potential risk to the user. Reprocessing, cleaning, disinfection and sterilization may compromise product characteristics which in turn create an additional risk of physical harm to or infection of the patient.

Adverse Events: The following events have been reported although their occurrence greatly depends on medical conditions of patient: infection, encrustation, obstruction, rupture, migration, bladder irritation symptoms, pain, hematuria, erosion.

See Instructions for Use for detailed information regarding warnings/precautions, adverse events prior to using this product. For further information contact Coloplast Corp at 1-800-258-3476 and/or consult the company website at www.coloplast.us

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

[Ostomy Care](#) / [Continence Care](#) / [Wound & Skin Care](#) / [Interventional Urology](#)

Coloplast Corp. Minneapolis, MN 55411 / Interventional Urology Surgical Support 1-800-258-3476

www.coloplast.us The Coloplast logo is a registered trademark of Coloplast A/S.

© 2019 Coloplast Corp. All rights reserved.

PM-08247 07.19