

Virtue[®]

Male Sling System

Providing *Elevation* and *Compression*

The Coloplast Virtue[®] Male Sling System is designed for stress urinary incontinence (SUI) in male patients. The Virtue[®] Sling is made of 100% monofilament polypropylene.

4-Arm Approach

- Designed to conform to patient's anatomy
- Ensures consistent placement with quadratic fixation
- Limits sling loosening and proximal or distal movement

Elevation and Compression

- Dual mechanisms of action contribute to patient continence³
- Transobturator arms provide proximal relocation through elevation of the bulbous urethra
- Pre-pubic arms provide urethral compression for increased resistance

Innovative Placement Technique

- Sub-urethral placement provides support to the bulbous urethra
- Reduced dissection leaves the urethra unexposed

89%

reduction
in pad weight

at 12 months for moderate range¹



92%

of patients
choose a male sling
over an AUS²



25%

of patients choose
a male sling vs. surgeon
AUS recommendation²



To learn more, visit ColoplastMD.com

 Coloplast

Ordering Information

Contact Customer Service at 800-258-3476

Product	Order Number
Virtue® Male Sling System with Alexis® Retractor	50020

Contains:

- 1 - Alexis® Retractor
- 1 - Virtue Sling
- 1 - Single-use Introducer



References

1. Comiter, Craig et al. The Virtue Sling—A New Quadratic Sling for Male Stress Incontinence—Results of a Multinational Clinical Trial. *Urology*, Vol. 84, Issue 1, Jul 2014.
2. Kumar, Angelish et al. Artificial Urinary Sphincter Versus Male Sling for Post-Prostatectomy Incontinence—What Do Patients Choose? *THE JOURNAL OF UROLOGY*, Vol. 181, 1231-1235, March 2009.
3. Comiter, Craig V. Virtue Quadratic Sling, a Product of Sling Evolution. AUA 2012.

Indications for Use: The Coloplast Virtue Male Sling System is an implantable, suburethral support sling indicated for the treatment of male stress urinary incontinence. Contraindications: The product is contraindicated for the following patients: with urinary tract infections or urinary tract obstruction; with blood coagulation disorders or prescribed anticoagulation therapy; with obstructive uropathy; under the age of 18.

Warnings: Carefully consider the risks and benefits of implantation in patients with the following: persistence of de novo incontinence; good bladder function and low volume detrusor overactivity; presence of bladder neck or urethral strictures; blood coagulation disorders, renal insufficiency due to urinary tract obstruction, or compromised immune systems. The patient should not have an active infection at the time of implantation.

Precautions: It is the responsibility of the surgeon to advise the prospective patients or their representatives, prior to surgery, of the possible warnings associated with the use of this product. This product should only be used by surgeons who are qualified to perform this type of surgery, and who are familiar with the use of non-absorbable mesh and the specific insertion technique for Virtue. A 6-month period of non-invasive treatment is recommended before a sling implant is considered for males with stress urinary incontinence.

Adverse Reactions: Local irritation, a foreign body response, wound dehiscence, extrusion, erosion, inflammation, fistula formation, urethral erosion, infection, bladder, urethra, vessel, nerve perforation. Other known risks of incontinence surgical procedures include sling migration, pain, transient or permanent retention, bladder outlet obstruction, and, continued stress urinary incontinence and persistent or new overactive bladder symptoms. The occurrence of these responses may require operative intervention with removal of part of or the complete sling.

See Instructions for Use for detailed information regarding warnings/precautions, adverse events prior to using this product. For further information contact Coloplast Corp at 1-800-258-3476 and/or consult the company website at www.coloplast.us.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Ostomy Care / Continence Care / Wound & Skin Care / Interventional Urology

Coloplast Corp. Minneapolis, MN 55411 / Interventional Urology Surgical Support 1-800-258-3476
www.coloplast.us The Coloplast logo is a registered trademark of Coloplast A/S.
© 2019 Coloplast Corp. All rights reserved.



PM-07790 04.19