

# Vortek<sup>®</sup>

## Double Loop Ureteral Stent

### *Essential value, reliable performance*

#### **The Value of Vortek<sup>®</sup>**

- Comprised of Vortek<sup>®</sup> material, designed for ease of insertion and compression resistance
- Ideal material for passing through tight ureters due to stenosis

#### **Dual-Layer Design Balances Insertion and Patient Comfort**

- Firm inner layer resists compression and allows for easy advancement into the ureter
- Thermosensitive outer layer softens at body temperature for increased patient comfort

#### **Smaller Size with Excellent Urinary Flow Rate**

- 4.8 Fr size has drainage capacity similar to a 6 Fr stent



## Ordering Information

Coloplast Urology Care Service Support 800-258-3476

French	Length	Order Number Without Guidewire	Order Number Seldinger Guidewire
4.8	12 cm	ACBM50	ACB1C0
	16 cm	ACBM51	ACB1C1
	20 cm	ACBM52	ACB1C2
	22 cm	ACBM57	ACB1C7
	24 cm	ACBM53	ACB1C3
	26 cm	ACBM54	ACB1C4
	28 cm	ACBM55	ACB1C5
	30 cm	ACBM56	
6	20 cm	ACBM62	ACB162
	22 cm	ACBM67	ACB167
	24 cm	ACBM63	ACB163
	26 cm	ACBM64	ACB164
	28 cm	ACBM65	ACB165
	30 cm	ACBM66	ACB166
7	20 cm	ACBM72	ACB172
	22 cm	ACBM77	ACB177
	24 cm	ACBM73	ACB173
	26 cm	ACBM74	ACB174
	28 cm	ACBM75	ACB175
	30 cm	ACBM76	ACB176
8	24 cm	ACBM83	ACB183
	26 cm	ACBM84	ACB184
	28 cm	ACBM85	ACB185
	30 cm	ACBM86	ACB186

**Indications for Use:** The double loop ureteral stents are used for: drainage of the upper urinary tract over fistulas or ureteral obstructions; healing of the ureter.

**Warnings:** Reuse of this single use product may create a potential risk to the user. Reprocessing, cleaning, disinfection and sterilization may compromise product characteristics which in turn create an additional risk of physical harm to or infection of the patient.

**Adverse Events:** The following events have been reported although their occurrence greatly depends on medical conditions of patient: infection, encrustation, obstruction, rupture, migration, bladder irritation symptoms, pain, hematuria, erosion.

See Instructions for Use for detailed information regarding warnings/precautions, adverse events prior to using this product. For further information contact Coloplast Corp at 1-800-258-3476 and/or consult the company website at [www.coloplast.us](http://www.coloplast.us).

**Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

[Ostomy Care](#)  
[Continence Care](#)  
[Wound & Skin Care](#)  
[Urology Care](#)

Coloplast Corp. Minneapolis, MN 55411 / Urology Care Surgical Support 1-800-258-3476  
[www.coloplast.us](http://www.coloplast.us) The Coloplast logo is a registered trademark of Coloplast A/S.

© 2018. Coloplast Corp. All rights reserved.

