

---

In response to your request for information, please review the following regarding Coloplast's saline filled testicular implants:

### **MR COMPATIBILITY**

The following testicular implants are fabricated with materials that are not conductive and are not magnetically reactive.

- P/N 450-1323
- P/N 450-1325
- P/N 450-1327
- P/N 450-1329

These implants are constructed of LSR 35 silicone rubber and fiberglass. These materials are non-conductive of electricity or magnetism.

Per ASTM F2503-08 items with non-conducting and non-magnetic construction may be deemed MR-Safe.


As the implants listed above meet these conditions Coloplast deems them to be MR-Safe.

### **LATEX CONTENT**

Coloplast's testicular prostheses are fabricated from materials that do not contain latex. Additionally, latex gloves are not worn by employees that manufacture these products.

For further information, please contact the Coloplast Surgical Support department at 800-258-3476.

Sincerely,



Joseph DuPay III  
Director of Global Research and Development  
Coloplast