


Dormia® DualPerc

PCNL Basket

Two solutions. One move.

Discover the power of simplicity with Dormia DualPerc.
Uniquely engineered to transform complex percutaneous
stone retrieval into a single, efficient action.


This is advanced care, *simplified.*



Maximize efficiency. Minimize waste.

Dormia DualPerc streamlines PCNL procedures with its 'one & done' approach

This single device embodies procedural efficiency and operational excellence, simplifying surgical processes. Its dual-purpose design minimizes waste and optimizes inventory.




Small steps. Big impact.

- **Reduced waste and inventory optimization:** Minimize waste and optimize inventory with Dormia DualPerc's dual-purpose design
- **Operational savings:** 2-in-1 device design switches to and from rigid to flexible basket configurations based on the surgeon's choice and needs during the procedure
- **Enhanced procedural versatility:** Smaller diameters than other percutaneous extractors allow Dormia DualPerc to adapt to a wider range of PCNL procedures

See it, remove it.

- **Mini-PCNL innovation:** The 4.8Fr Dormia DualPerc is the first stone retrieval basket designed for mini-PCNL
- **Enabling stone-free rates:** Dual-purpose rigid and flexible design enables immediate removal of residual fragments



It's not just efficient — it's your practice, **optimized.**

Together in care_

At Coloplast, we're dedicated to improving patient care with innovative solutions that simplify complex procedures. Our commitment goes beyond products; Together, we're achieving healthcare excellence.



Color coded handle
to easily differentiate diameter size

Ordering Information

Coloplast Interventional Urology Service Support 800-258-3476

Item	Description	Rigid Shaft Length	Rigid Shaft Diameter	Flexible Shaft Length	Flexible Shaft Diameter	Basket Size	Sales UOM	EA/Sales UOM
EXP448	Dormia DualPerc Stone Extractor 4.8Fr	35 cm	4.8Fr	64 cm	3Fr	12.5 mm	BX	5
EXP470	Dormia DualPerc Stone Extractor 7Fr	35 cm	7Fr	64 cm	3Fr	12.5 mm	BX	5

DORMIA® DUALPERC BRIEF STATEMENT

Intended Use

Dormia DualPerc is intended to be used for extraction of stones in the kidney and bladder.

Intended Medical Indication

Dormia DualPerc is an interventional endourological device used to remove stones and stone debris following lithotripsy, located in kidney (calix, pelvis, ureteropelvic junction...) through antegrade percutaneous access or in the bladder through retrograde natural access.

Contraindications

Contraindications to the medical device

The device is not intended for use other than for endourological procedure. Do not use, when, in the judgement of the physician, such a procedure would be contrary to the best interest of the patient. The device contains nitinol materials. Do not use in case of patient's known allergy to device materials.

Contraindications to the endourological procedure

Untreated urinary tract infection. Uncontrolled haemostasis disorder. Malignant kidney or urinary tract tumor. Percutaneous NephroLithotomy (PCNL) procedure is not indicated in:

- Pregnancy unless exceptional situation
- Bleeding disorders

Warnings and Precautions

This device must only be used by trained and experienced physicians. Physicians must inform patients of the possible undesirable side effects.

Potential Complications

Adverse events involving patients have been reported, which may be related to the procedures used or problems occurring with the device (such as detachment of a part of the device, breakage, or inability to withdraw the basket), especially if the operating procedure and the warnings above are not observed, although their occurrence greatly depends on patients' medical conditions. Adverse events include but are not limited to: Tissue trauma (mucosal abrasion, perforation), infection (e.g. urinary tract infection, pyelonephritis, severe infection...), pain, bleeding, haematuria, foreign object in body, device embedded in tissue, prolonged procedure and stone impaction.

Advice to the Patient

The patients should be educated about their interventional procedure. The patients should be advised in terms that they understand to inform the attending physician immediately if any abnormality or adverse event occurs (as haematuria, infection, pain).

The risks and benefits of using Dormia® DualPerc should be considered in patients.

The information provided is not comprehensive with regard to product risks. For a comprehensive listing of indications, contraindications, warnings, precautions, and adverse events refer to the product's Instructions for Use. Alternatively, you may contact a Coloplast representative at 1-800-258-3476 and/or visit the company website at www.coloplast.com.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Minneapolis, MN PM-28278